COVID-19 & Influenza A+B Antigen Combo Rapid Test For Self-Testing Package Insert

CE 2934 IVD REF FCO-6032H English

INTENDED USE

The COVID-19 Influenza A+B Antigen Combo Rapid Test is a single-use test kit intended for qualitative detection of nucleocapsid protein antigens of Influenza A and B viral antigens and COVID-19 Antigen from nasal swabs in individuals aged 12 and older. This test is intended for home use with self-collected nasal swab samples in individuals aged 12 and older. Sampling from anyone under the age of 12 should be performed only under the supervision and guidance and assistance of their guardian. People who are unable to carry out the test on their own should seek support.

The test is intended for asymptomatic individuals within 7 onset days or asymptomatic individuals in contact with persons who have been diagnosed positive or with individuals suspected to be infected. Positive results are indicative of the presence of influenza and SARS-CoV-2. Individuals who are tested positive should continue to wear additional personal protective equipment. Positive results do not rule out bacterial infection or co-infection with other viruses. Negative results do not preclude influenza and SARS-CoV-2 infection. Individuals who test negative and continue to experience influenza or COVID-19 symptoms should seek follow-up care from their healthcare provider.

PRINCIPLE

The COVID-19 Influenza A+B Antigen Combo Rapid Test is an immunochromatographic membrane assay that uses highly sensitive antibodies to detect SARS-CoV-2, Influenza A and B nucleocapsid protein from nasal swab samples. Influenza A and B specific antibodies are immobilized onto the test region of the membrane and combined with other reagents to conduct a test strip. The test is designed to detect nucleocapsid protein antigens in nasal swab specimens, which is different from the mutation sites have occurred in the spike protein, so it is theoretically able to detect variants including those in UK, India, South Africa and Brazil.

KIT COMPONENTS

Component	1 Test Kit	2 Tests Kit	5 Tests Kit	10 Tests Kit	20 Tests Kit	25 Tests Kit
COVID-19 Influenza A+B Antigen Combo Test	1	2	5	10	20	25
Extraction tube with buffer	1	2	5	10	20	25
Sterilized nasal swab	1	2	5	10	20	25
Workstation	/	/	/	/	/	/
Package insert	1	1	1	1	1	1

ADDITIONAL SPECIAL EQUIPMENT

- Timer
- If the test result is positive
 - Contact your doctor / general practitioner or the local health department immediately
 - Comply with local guidelines for self-isolation
 - To have a PCR confirmatory test performed
- Warnings and Precautions
 - Do not use after expiration date. Do not use if kit is damaged or open. Do not reuse the tests.
 - Do not eat, drink or smoke in the area where the specimens or kits are being handled.
 - Handle all specimens as if they contain infectious agents. Discard the using testing materials in accordance with local regulations.
 - The extraction buffer contains a salt solution if the solution contacts the skin or eyes, flush with clean water. Avoid contact with the eyes. If contact with the eyes occurs, flush thoroughly with water and give plenty of water to dilute the substance. If any discomfort, seek medical attention immediately.
 - Children and elderly please use the test under the guidance.

STORAGE AND STABILITY

Store unused test devices unopened at 4°C-30°C. If stored at 4°C-8°C, ensure that the test device is brought to room temperature before use. The test device is stable through the expiration date under the conditions. Do not freeze the past the kit's expiration date.


TEST PROCEDURE

- Open the kit box. Check the components before use. Please read all instructions carefully before you begin.
- Prepare before sampling:
 - Get a flat surface ready, like a table. Make sure it is clear, clean and dry.
 - Take off hand jewelry.
 - Wash your hands for 20 seconds. Use soap and water, or hand sanitizer. Dry your hands using clean disposable paper towels.
 - For best protection and avoid cross-contamination, disposable gloves, masks and eye protection (not provided in the package) are recommended.
- Cross-reactivity: Cross-reactivity studies are performed to demonstrate that the test does not react with the following microorganisms in the table below at concentration of 1x10⁷ TCID₅₀ ml for viruses and 1x10⁷ CFU/ml for bacteria.

Substrate	Influenza A/B			COVID-19		
	Positive	Negative	Total	Positive	Negative	Total
Salivary	42	42	84	95	95	190
Urine	42	42	84	95	95	190
Stool	42	568	610	110	455	565
Relative Specificity	100.00% (95%CI:91.59% - 100.00%)	86.38% (95%CI:78.51% - 92.16%)		100.00% (95%CI:99.35% - 100.00%)	100.00% (95%CI: 98.19% - 100.00%)	
Overall Agreement	100.00% (95%CI:99.35% - 100.00%)	100.00% (95%CI: 98.19% - 100.00%)		100.00% (95%CI: 99.35% - 100.00%)	100.00% (95%CI: 98.19% - 100.00%)	

- Cross-reactivity: Cross-reactivity studies are performed to demonstrate that the test does not react with the following microorganisms in the table below at concentration of 1x10⁷ TCID₅₀ ml for viruses and 1x10⁷ CFU/ml for bacteria.

- Open the sealed pouch and take out the test cassettes. For best results, the test should be performed in one hour.
- Insert the tube vertically up to the sample well.
- Add 3 drops specimen to the each sample well by gently squeezing the sides of the tube, then start the timer.

INDEX OF SYMBOLS

Symbol	Description	Symbol	Description
🚫	Do not reuse	📄	For in vitro diagnostic use only
🌡️	Store between 4-30°C	📄	Consult instruction for use
☀️	Keep away from sunlight	📄	Lot number
🕒	Use-by date	📄	Contents sufficient for <=> tests
🏭	Manufacturer	📄	Date of Manufacture
📄	Sterilized using ethylene oxide	📄	Sterilized using irradiation
📄	CE Mark	📄	CE Mark
📄	Authorized Representative in the European Community	📄	Contains biological material of animal origin

Result	Interpretation
Two distinct colored lines appear in the left window. One colored line should be in the control region (C) and another colored line should be in the Influenza A region (A).	Positive
Two distinct colored lines appear in the left window. One colored line should be in the control region (C) and another colored line should be in the Influenza B region (B).	Positive
Two distinct colored lines appear in the left window. One colored line should be in the control region (C) and another colored line should be in the Influenza A region (A) and Influenza B region (B).	Positive
No colored line appears in the control region (C) of the left window. No apparent colored line appears in the test region (B/A).	Negative
No colored line appears in the control region (C) of the left window. No apparent colored line appears in the test region (B/A).	Negative
No colored line appears in the control region (C) of the left window. No apparent colored line appears in the test region (B/A).	Negative

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TEST PROCEDURE

- Open the kit box. Check the components before use. Please read all instructions carefully before you begin.
- Prepare before sampling:
 - Get a flat surface ready, like a table. Make sure it is clear, clean and dry.
 - Take off hand jewelry.
 - Wash your hands for 20 seconds. Use soap and water, or hand sanitizer. Dry your hands using clean disposable paper towels.
 - For best protection and avoid cross-contamination, disposable gloves, masks and eye protection (not provided in the package) are recommended.
- Cross-reactivity: Cross-reactivity studies are performed to demonstrate that the test does not react with the following microorganisms in the table below at concentration of 1x10⁷ TCID